Cleaning & Disinfection

Adequate cleaning is considered the most essential step in the reprocessing cycle of sterile goods.

Preface

The Instrument Reprocessing Working Group (AKI) was founded in 1976 with the aim of developing and publishing expertise and practical information on the safe reprocessing of instruments in healthcare with the aim of retaining their value. Having been established for four years as of 2016, AKI has not only become part of the history of "instrument reprocessing", but also embodies the present, serving as an inspiration for the future of our speciality field as an example of multidisciplinary and interactive cooperation.

Efforts to ensure the independence of the departments for sterile supplies over these past decades, now known as the "Reprocessing Unit for Medical Products (RUMP)", were an important step towards establishing a scientific basis and much improved reprocessing. This was later followed by the industrialization of the department (i.e. viewing the sterile supplies department as an industrial production unit), assisted by the introduction of quality management systems and IT systems.

Today, we have reached the next phase, defined by a more unified perspective of "reprocessing". A sterile medical product is the sum of a series of actions, of which the sterilization process itself is only a part. Sterility is the result of integrity across processes, which also requires "all-inclusive" responsibility to be assumed by the reprocessing department for this overall process. This responsibility no longer ends at the exit door of the Reprocessing Unit for Medical Products – today, it reaches as far as the operating theater, the patient's bed, and the treatment room of the clinic.

AKI's great achievement lies in always having remained within its essence. Its various publications have made a significant contribution to establishing a robust basis for the reprocessing of healthcare instruments, offering practical solutions to everyday problems faced by sterile supplies departments or dental practice. For instance, these may be problems relating to water reprocessing and water quality, to the selection and dosage of cleaning agents, to the structure of programs of various devices, to the quality of the instruments themselves, and also to the organization of work flows, which has a major impact on the quality of the end product and on its long-term value. These are recurring problems for which even the latest apps cannot help to find definitive solutions. Ultimately, they require daily monitoring of the various parameters and timely intervention in the event of significant discrepancies. They require the ongoing vigilance of the staff of the Reprocessing Unit for Medical Products.

The influence that AKI's publications have had and continue to have on modern sterilization is impossible to overestimate. AKI since offers

www.a-k-i.org

Instrument Reprocessing

Reprocessing
of Instruments to Retain Value

Working Group
Instrument Reprocessing

11

5. Preparation for Cleaning and Disinfecting

The first steps in a proper reprocessing cycle are taken in the operating theater. Coarse contaminants, residues from hemostatics, skin disinfectants and lubricants, as well as caustic drugs should, wherever possible, be removed before the instruments are set down.

Never immerse stainless steel instruments in an isotonic solution (such as physiological saline solution). This is because prolonged instrument contact with saline solution leads to pitting and stress corrosion cracking.

Waste, skin disinfectant residues, saline solutions etc., may not be put in disposal containers.

In hospitals with a Reprocessing Unit for Medical Products, self-contained systems are used to transport contaminated medical products from the operating theaters and wards to the unit. Dry disposal always preferable as a more conservative method aimed at retaining value.

Wet disposal involves damp cloths being placed on the instrument trays during transport or the use of an appropriate foam spray. The damp conditions increase the risk of corrosion and may encourage the growth of microorganisms.


Because of the corrosion risk and the cleaning factors, long intervals between instrument use and reprocessing (e.g. overnight or over the weekend) should be avoided. Field experience has shown that in the case of dry disposal, intervals of up to 6 hours pose no problem. Contamination and pre-cleaning are factors that have a critical impact here.

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Water treatment

	Drinking water		Softened water		Demineralized water
Evaporation residue (ppm)	500		530		5
Elec. conductivity ($\mu\text{S}/\text{cm}$)	650		700		3
Total hardness ($^{\circ}\text{d}$)	14		< 0.1		< 0.1
Sodium salts (mg/l)	20		160		< 1
Chlorides (mg/l)	40		40		< 1
Silicates (ppm SiO_2)	12		12		< 0.1
pH value	6.7		8		5.5

***Chloride content do not exceed 100 mg/l at room temperature.**

Spotting caused by Lime


Stains on instruments. Cause:
Dried water with high salt content

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

Due to the fact that silicates do not increase the conductivity of the water, the presence of silicates may be overlooked easily.

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Spotting caused by silicon dioxide in the steam condensate.

Chloride pitting

Chlorides are dangerous

Chloride-induced pitting on instrument.

Chlorides

Chlorides dissolved in the water are particularly critical substances, as they tend to cause pitting even on stainless steel instruments if present in higher concentrations.

The danger of chloride-induced pitting generally rises with:

- An increase in the chloride content,
- An increase in temperature,
- Declining pH value,
- Increasing exposure time,
- Insufficient drying,
- Concentration of chloride resulting from adherence of dry residues to instrument surfaces after evaporation.

While the causal relationships between the chloride content of the water and pitting are not always predictable, laboratory testing has shown sign of corrosion on instruments after just two hours with a chloride content of 100 mg/l at room temperature. As chloride concentrations increase, the risk of pitting also increases rapidly.

Chlorides are dangerous

Scanning electron microscope image, chloride-induced pitting

5. Preparation for Cleaning and Disinfecting

DR. WEIGERT

The first steps in a proper reprocessing cycle are taken in the operating theater. Coarse contaminants, residues from hemostatics, skin disinfectants and lubricants, as well as caustic drugs should, wherever possible, be removed before the instruments are set down.

Never immerse stainless steel instruments in an isotonic solution (such as physiological saline solution). This is because prolonged instrument contact with saline solution leads to pitting and stress corrosion cracking.

Careless dropping can also damage instruments. For example, the hardened (tungsten carbide) tips of scissors may come off, or small clamps may be bent. To avoid damage, always put your instruments down carefully after use. Do not overload instrument trays.

Waste, skin disinfectant residues, saline solutions etc., may not be put in disposal containers.

In hospitals with a Reprocessing Unit for Medical Products, self-contained systems are used to transport contaminated medical products from the operating theaters and wards to the unit. Dry disposal always preferable as a more conservative method aimed at retaining value.

Wet disposal involves damp cloths being placed on the instrument trays during transport or the use of an appropriate foam spray. The damp conditions increase the risk of corrosion and may encourage the growth of microorganisms.

Because of the corrosion risk and the cleaning factors, long intervals between instrument use and reprocessing (e.g. overnight or over the weekend) should be avoided. Field experience has shown that in the case of dry disposal, intervals of up to 6 hours pose no problem. Contamination and pre-cleaning are factors that have a critical impact here.

Corrosion caused by immersion in physiological salt solution over a period of several hours

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Disinfecting cleaning

When performing disinfecting cleaning, please also note that the disinfecting capability should have been proven under "dirty conditions" (high protein load) in accordance with European (EN) standards or corresponding national regulations. When using disinfecting cleaning agents, the producer's specifications regarding concentration, temperature, exposure time and material compatibility must be observed. The cleaning/disinfecting solutions used should be freshly prepared on a daily basis. If it becomes visibly contaminated, it is advisable to prepare fresh solutions at even shorter intervals.

Under dirty conditions (non-pre-cleaned instruments) 20 °C	
bactericidal, yeasticidal (VAH)	10 ml/l (1.0 %), 15 min
bactericidal (EN 13727, EN 14561), yeasticidal (EN 13624, EN 14562)	10 ml/l (1.0 %), 5 min or 5 ml/l (0.5 %), 15 min
Active against enveloped viruses (incl. HIV, HBV, HCV) (EN 14476, EN 17111)	10 ml/l (1.0 %), 5 min or 5 ml/l (0.5 %), 15 min
Active against enveloped viruses (incl. HIV, HBV, HCV) (DVV/RKI)	15 ml/l (1.5 %), 5 min or 10 ml/l (1.0 %), 15 min
Rotavirus (EN 14476, EN 17111)	10 ml/l (1.0 %), 15 min
Polyoma SV 40 (EN 14476, EN 17111)	20 ml/l (2.0 %), 15 min or 15 ml/l (1.5 %), 30 min

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

Chloride sources in the instrument usage and processing cycle:

- Fresh-water chloride content (depending on the source of the supply)
- Insufficient demineralization of the water used for the final rinse and steam sterilization
- Regeneration salt carryover, leakage or spillage from ion exchangers used for water softening
- Use of agents not permitted for or incorrectly used during reprocessing
- Isotonic solutions (e.g. a physiological salt solution), etchants and drug residues
- Organic residues (body fluids such as blood, chloride content 3200-3550 mg/ltr, saliva, sweat) dried on the surfaces

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


Neodisher Prestop

neodisher PreStop – new Corrosion – inhibition test!

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
Safety data sheet in accordance with regulation (EC) No 1907/2006

neodisher MultiZym

Version: 2 / GB

Replaces Version: 1 /
GB

Date revised: 11.03.2015

Print date: 01.04.15

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

neodisher MultiZym

1.2. Relevant identified uses of the substance or mixture and uses advised against

Use of the substance/preparation

Washing and cleaning products (including solvent based products)

1.3. Details of the supplier of the safety data sheet

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(Toxicology Test EN ISO 10993-1 biological assesment of medical products)

Toxicological Profile for Neodisher LM 2	page 2 of 4
---	-------------

Genotoxicity (Mutagenicity)

The product does not contain any ingredients classified or known for such properties(4,8).

Reproductive toxicity

The product does not contain any ingredients classified or known for such properties (4,8).

Carcinogenicity

The product does not contain any ingredients classified or known for such properties(4,8).

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- Ultimate or inherent degradation

OECD 301 A-F:

All organic components are readily biodegradable.

The surfactants contained in the product are biodegradable according to the requirements of the Detergent Directive 648/2004/EC (10).

Bioaccumulation

The product does not contain any bioaccumulating substances.

German water hazard class 1

Slightly harmful to aquatic organisms.

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